
READ ALL OF THIS LEAFLET CAREFULLY BEFORE YOU START USING THIS SOLUTION

Keep this leaflet. You may need to read it again.


INSTRUCTION FOR USE / HANDLING:

The product is presented in a two-compartment bag separated by a peel seal. Aseptic technique should be used throughout the administration to the patient. Use only if the solution is clear, all seals are intact and the overwrap is undamaged. All seals, including the peel seal between the two compartments must be intact. If leakage is discovered, discard the solution immediately since sterility can no longer be assured. Not for direct infusion.

MIXING


I Immediately before use, remove the overwrap from the bag and mix the solutions in the two different compartments. Hold the small compartment with both hands and squeeze it until an opening is created in the peel seal.

II Squeeze with both hands on the large compartment until the peel seal between the two compartments is entirely open. Shake gently to mix. The solution is now ready for use and can be hung on the equipment.

CONNECTING

The dialysate line may be connected to bag through either the luer or the spike connector.

III The luer connector is a needle-less and swabbable connector. Remove the cap with a twist and pull motion, and connect the male luer lock on the dialysate line to the female luer receptor on the bag. Ensure that the connection is fully seated and tighten. The connector is now open. Verify that the fluid is flowing freely during use. When the dialysate line is disconnected from the luer connector, the connector will close and the flow of the solution will stop.

If the **spike connector** is used, first remove the snap-off cap. Then introduce the dialysate line spike through the rubber septum of the bag connector. Ensure that the spike is fully inserted and verify that the fluid is flowing freely.

The reconstituted solution should be used up within 24 hours of removing the overwrap. The reconstituted solution is for single use only. Discard any unused solution immediately after use.

During treatment, hemodynamic status, fluid balance, electrolyte and acid-base balance should be closely monitored throughout the procedure.

WHAT DOES PRISMASATE SOLUTION CONTAIN?

The product is presented in a two compartment bag containing in the smaller compartment A, an electrolyte solution, and in the larger compartment B, a buffer solution. Solution formulations for each compartment are listed on the bag. The final reconstituted solution is obtained after breaking the peel seal between the two compartments and mixing both solutions. Refer to bag labeling for composition of the final reconstituted solution.

PACKAGING:

2 bags at 5000 mL each per case.

WHAT IS PRISMASATE SOLUTION AND WHAT IS IT USED FOR?

The product is a sterile dialysis solution intended for treatment of acute kidney disease (renal failure) using Continuous Renal Replacement Therapies, such as continuous hemodialysis and hemodiafil-

tration aimed at normalizing the composition of the blood.

The product is a sterile solution packaged in a two-compartment bag. The electrolyte solution (in the small compartment A) must be mixed with the buffer solution (in the large compartment B) before use to obtain the final solution suitable for the treatment.

The product may also be used in case of drug poisoning with dialyzable or filterable substances.

The product should only be used by or under the prescription and direction of a physician competent in intensive care treatment using Continuous Renal Replacement Therapy.

Federal Law (USA) restricts this device to sale by or on the order of a physician.

WHAT SIDE EFFECTS CAN PRISMASATE SOLUTION CAUSE?

When continuous hemodialysis and hemodiafiltration are performed correctly, side effects are uncommon.

Some side effects may occur, including nausea, vomiting, muscle cramps and low blood pressure (hypotension).

STORAGE AND EXPIRY DATE

Do not store below +4°C/39°F. Do not use after the expiry date shown on the label and the packaging.

From a chemical point of view, once the overwrap is removed, the product should be used immediately. Stability of the product after removal of the overwrap is 24 hours including the time of treatment. Other in-use storage times and conditions are the responsibility of the user.

This product is not made with natural rubber latex.

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